

# Guardian® Class I Biological Safety Cabinets



MSC 1200 C1

**BRITISH  
MANUFACTURED**



**FIVE YEAR  
WARRANTY**



**MAXIMUM  
OPERATOR  
SAFETY**



**GUARANTEED  
PROTECTION**



**The Guardian® Class I Biological Safety Cabinet offers operator and environmental protection from biological materials and activities involving animal pathogens and genetically modified organisms.**

Designed to incorporate the latest features to BS EN 12469:2000 Standards, the range of Microbiological Safety Cabinets offer total practicality and suitability for a huge variety of applications.

The units are provided as a re-circulating, or optionally ducted, system and utilise HEPA filtration to offer environmental and operator protection. The H14 HEPA's are rated at 99.997% efficient at 0.3µm particulate, the filter seals under negative pressure to prevent leaks into the cabinet.

The MSC1200C1 can be fitted with Carbon Filters to allow use with a multitude of chemical compounds, including formalin. Universal ports allow connection to proprietary VHP and formalin systems to ensure safe and complete disinfection.

## Size Specification (W x D x H)

**MSC 1200 C1**

1200mm x 700mm x 1350mm \*

\*Add 15mm height for Outlet Carbon Filter, if fitted.